


## Formaldehyde Neutralizing Pads

PART #	SIZE	PACKAGING
FNP-0810B-400	8" x 10" Blue	400 pads per case
FNP-1620B-100	16" x 20" Blue	100 pads per case
FNP-1212-200	12" x 12" White Wipe	200 wipes per case



FNP-0810B-400 & FNP-1620B-100



FNP-1212-200

CLINICAL INNOVATION • SUSTAINABLE TECHNOLOGY • ULTIMATE PERFORMANCE

### Product Details—Neutralize Dangerous Formaldehyde Spills

- **Treated Soaker Fabric**—Absorbs, contains, and neutralizes formalin spills and waste; and chemically neutralizes the formaldehyde.
- **Workplace Safety**—Provides simple means for formalin spill containment and disposal
- **Neutralize**—Chemically neutralizes both formaldehyde and glutaraldehyde
- **Hydroscopic Design**—Helps absorb harmful vapors from the air (white wipe)

### Manage Risk and Comply with OSHA Standards

- **Reduces Exposure**—Exposure to formaldehyde; formaldehyde is a toxin and carcinogen
- **Minimize Contact With Spill And Vapors**—Which can cause skin, eye, and respiratory irritation

### Eco-Friendly—Minimizing Chemicals of Concern and Waste

- PVC-free and latex-free
- Safe for incineration