


| PART # | SIZE | PACKAGING |
|----------|-----------|------------------|
| ESP-3630 | 36" x 30" | 25 pads per case |
| ESP-3640 | 36" x 40" | 25 pads per case |
| ESP-3072 | 30" x 72" | 25 pads per case |
| ESP-4672 | 46" x 72" | 10 pads per case |



Multiple Applications

- *Fluid Intense Procedures
- *In Front of Scrub Sinks
- *Organ and Tissue Recovery

Product Details—Outstanding Absorption, Outstanding Value

- **Absorbent**—Absorbs approximately 1 gallon of fluid per square yard
- **Durable**—Strong, quilted design resists pilling, holds up to foot scuffing and stool wheels
- **Stable**—Non-skid barrier keeps the floor under the pad clean and dry for easier cleanup and improved safety

Manage Risk and Comply with OSHA Standards

- **Compliance**—Aides compliance to the OSHA Bloodborne Pathogen Standard
- **Reduce Exposure**—To fluids per OSHA Standards for employee safety
- **Avoid Trips, Slips and Falls**—Improving safety of healthcare professionals

Eco-Friendly—Minimizing Chemicals of Concern and Waste

- PVC-free and latex-free
- Manufactured with recycled materials
- Low VOC emissions during incineration, for crematorium compatibility