

LDI medical™

Sustainable Technology · Clinical Innovation · Ultimate Performance

Positioners for positive outcomes ...

EnviroSoft® Reusable Surgical Positioning Devices: New patent-pending high-performance tools can be safely cleaned + used up to 50 times with confidence

Protect your patients from risks associated with decubitus ulcers and nerve damage during surgery while reducing waste, thanks to a new patent-pending product line engineered by the LDI Medical team.

EnviroSoft® Reusable Surgical Positioning Devices feature a high-performance, strategic foam that's 100% impervious to fluid and easily cleanable with standard surgical wipes, allowing safe use up to 50 times before disposal. **Discover your new choice for reusable positioning devices.** Available in 3 configurations for comfort, stability, and risk management. Conform to your protocols that are required for safe pressure reduction, weight distribution and positioning – without the hassles or costs associated with super-expensive gel positioners or wasteful single-use disposables.

The Science Behind the Solution

- ✓ **Advanced Technology + Materials:** Durable, unique foam material maximizes support of soft tissue and minimizes pressure around bony prominences, reducing risk of pressure ulcers. Patent-pending **EnviroSoft® Surgical Positioning Devices** are 100% waterproof (impervious to fluid to the core), even when nicked or cut.
- ✓ **Cost-Effective, Cleanable for Long-Lasting Performance:** Unique supportive foam construction and thoughtful seamless design won't "bottom out" like single-use positioners and eliminates the concerns (and costly repairs) of expensive gel and reusable foam positioners that leak and harbor bacteria when nicked or cut. **Quickly clean EnviroSoft® positioners with standard surgical wipes.**
- ✓ **Eco-Friendly:** 100% PVC-free and latex-free. No flame retardants or antimicrobial coatings for allergen reduction. Safely reusable (up to 50 times with proper care) reduces waste.

EnviroSoft® Cleanable, Reusable Surgical Positioning Devices include Head Donut, Chest Roll and Arm Cradle.

LDI Medical - LDI Solutions

Contact Jack O'Callaghan • (866) 332-0700
www.ldimedical.com • sales@ldisolutions.com

