

EnviroSoft® Surgical Positioning Devices FAQ's

- **What is different about the LDI EnviroSoft® Positioners?**
 - Reusable - EnviroSoft® Surgical Positioners can be used up to 50 times.
 - Easy to Clean - Clean with standard surgical cleaning wipes. A list of approved/tested cleaners is available.
 - Impervious to fluid - Unlike most reusable positioners that can leak or absorb liquid when nicked or cut, the EnviroSoft® Surgical Positioners can still be used and cleaned if nicked or cut.
- **What are they made of?**
 - EnviroSoft® Surgical Positioners are made of a unique foam that does not absorb fluid, it is completely impervious to fluid.
- **Are they fluid-proof?**
 - Yes.
- **What types of positioners are available?**
 - Supine head positioner
 - Arm cradle
 - Chest roll
 - Prone head positioner
 - Heel cup
- **What can be used to clean the positioners?**
 - Standard surgical cleaning wipes.
 - A list of tested/approved cleaners and disinfectants is available. If there is a cleaner that you use and we have not tested it, we will test it for you.
- **How many times can they be used?**
 - Reusable up to 50 times.