

## EnviroSorb® Absorbent Skin Reconstruction Pack FAQ's

- **What is the Skin Reconstruction Pack used for?**
  - The SRK is designed to be used post recovery to absorb fluids and reconstruct the torso, buttock and shoulders. It can be customized but cutting away excess material to fit any donor and the absorbent materials help manage fluid during transport to the funeral home.
- **What is the material made of?**
  - High density polyethylene, PVC-free and latex free, Made with over 80% post-consumer recycled material
- **How much fluid can the SRP hold?**
  - The SRP can absorb one gallon of fluid per square yard. That means each SRP can absorb up to eight gallons of fluid. That will accommodate most donors.
- **Why is the collar around the opening for the head thicker?**
  - The additional, super-absorbent material is designed to stop any fluid ingress due to gravity. It will help keep fluid from staining the face during transport.
- **What is the pad in the middle designed for?**
  - The super absorbent pad in middle is designed to collect increased fluids in this location.
- **How many per case?**
  - 10 Skin Reconstruction Packs per case.
- **Some donors are smaller, this seems to have too much material, what should I do with that?**
  - Skin Reconstruction Packs are made to cut unwanted material and relocate to areas if needed.
- **Are they leak proof?**
  - The Skin Reconstruction Pack is leak resistant. For maximum leak protection use with The LDI Mort-Port™ Unionall.
- **Are there different sizes?**
  - Not currently.